

INDICATIONS

NuStat is a single-use hemostatic wound dressing applied externally with mechanical compression to temporarily control bleeding in lacerations, punctures, abrasions, surgical wounds (operative, postoperative, dermatological, etc.) and traumatic injuries.

DEVICE DESCRIPTION

NuStat Flex[®] XR is a sterile hemostatic temporary surgical dressing comprised of a knitted blend of silica and cellulose based fibers, and an X-ray detectable element.



CONTRAINDICATIONS

Do not leave NuStat in place for more than 24 hours

PRODUCT CODE	SIZE	UNITS PER BOX
XR-3000097333	12"x12"	10
XR-3000097317	3"x48"	10
XR-3000097320	4"x48"	10
XR-3000097331	6"x60"	10

WARNINGS

NuStat Flex[®] XR is not intended as a substitute for meticulous surgical technique, the use of sterile procedures, & the proper application of ligatures or other conventional procedures to control bleeding.



See IFU



Do not use if package is damaged



Keep dry



Do not reuse



Sterilized using radiation



Do not resterilize



Caution: U.S. Federal law restricts this device to sale by or on the order of a physician



Not made with natural rubber latex

INSTRUCTIONS FOR USE

This Instruction for Use will be used in conjunction with NuStat Flex® XR. See OVIK Health product codes above.

1. As desired, prewet NuStat Flex® XR with sterile saline/water.
2. Apply NuStat Flex® XR directly to the inner surface of the wound.
3. Apply firm pressure so NuStat Flex® XR makes direct contact with the bleeding area.
4. Maintain pressure and secure bandage if necessary.
5. As desired, rewet while in use to ease removal.
6. If bleeding continues, apply new NuStat Flex® XR as directed above.

PRECAUTIONS

- NuStat Flex® XR dressing is supplied as a sterile, single use dressing. Do not re-use. Do not re-sterilize.
- If bleeding resume after careful removal of the dressing, the use of standard surgical techniques may be required to adequately address the source of bleeding.

Ovik HEALTH

Spartanburg, SC USA

800.431.6686

ovikhealth.com

 Developed and manufactured by:

BEEKEN
BIOMEDICAL®

Newton, MA 02459 USA

888.497.7376



STERILE R

