
TWO-LAYER COMPRESSION KIT WITH ZINC

Instructions For Use

LAYER 1
← Absorptive
Foam Layer
with Zinc

LAYER 2
→ Cohesive
Compression Layer

- | | |
|--|--|
| (ENG) Two-Layer Compression with Zinc | (MS) Mampatan Dua Lapisan dengan Zink |
| (AR) ضغط ثنائي الطبقات مع الزنك | (NL) Tweelaagse compressie met zink |
| (CS) Dvouvrstvá komprese se zinkem | (NO) To-lags kompresjon med sink |
| (DA) To-lags kompres med zink | (PL) Dwuwarstwowy system uciskowy z cynkiem |
| (DE) 2-Schicht-Komprimierung mit Zink | (PT) Compressão de duas camadas com zinco |
| (ES) Compresión de dos capas con zinc | (RU) Двухслойная компрессия с цинком |
| (FI) 2-kerroksinen kompressiosidos sinkillä | (SK) Dvojvrstvová kompresia so zinkom |
| (FR) Compression double couche avec zinc | (SL) Dvoslojna kompresija s cinkom |
| (HE) קומפרסיה דו-שכבתית עם אבץ | (SV) Tvålagers komprimering med zink |
| (HR) Dvoslojna kompresa s cinkom | (ZH-CN) 锌两层压缩 |
| (IT) Compressione a due strati con zinco | (ZH-TW) 兩層壓縮含鋅 |

Description

These instructions are for both CoFlex TLC Zinc and CoFlex TLC Zinc Lite.

CoFlex TLC Zinc (*Standard & Lite*) is a two-step non-latex compression system. Layer 1 is a foam roll with Zinc that provides comfort and absorbs drainage. Layer 2 is a short stretch, moisture-resistant compression bandage that provides the proper levels of sustained working pressures to promote venous return. This system is designed to work together and provide continuous restorative compression for up to one week. A nylon stocking is included for ease of use and patient comfort.

Indications for Use

CoFlex TLC Zinc (*Standard & Lite*) is designed to treat venous leg ulcers and other related conditions.

CoFlex TLC Zinc (*Standard*) can be used for patients with an ABI greater than or equal to 0.8.

CoFlex TLC Zinc (*Lite*) can be used for patients with an ABI greater than or equal to 0.5.

Contraindications

CoFlex TLC Zinc (*Standard*): ABI less than 0.8.

CoFlex TLC Zinc (*Lite*): ABI less than 0.5.

Precautions & Warnings

- Should be used under the care of a licensed healthcare professional.
- Measure the Ankle Brachial Index prior to application.
- CoFlex TLC Zinc (*Standard & Lite*) Layer 1 is not a primary wound dressing. Follow appropriate protocol for the wound condition before applying.
- CoFlex TLC Zinc (*Standard & Lite*) Layer 2 may initially be applied with less tension allowing the patient to adjust to or tolerate the short stretch compression system.
- At the discretion of a healthcare professional, patients or their care providers can be trained to apply the bandage.
- Instruct patient to contact physician or nurse if they experience pain, numbness, tingling, discoloration or if swelling occurs.
- Product is not made with natural rubber latex.
- Single use item. If reused, intended restorative compression cannot be assured and increases potential for infection.
- Layer 1 of this product contains aloe.
- If redness or irritation occurs, discontinue use.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


Directions for Use

General Instructions:

Position the foot at a 90° angle.

Layer 1 - Foam Layer:

- (1) Begin wrapping Layer 1 on the top of the foot, at the base of the fifth digit in a spiral technique (figure 1).
- (2) Continue spiral wrapping across the top of the foot then proceed to wrap 1 figure of 8 around the back of the heel/ankle. The plantar surface of the heel may be partially exposed (figure 2).
- (3) Continue spiral wrapping up the leg. Apply with a 50% overlap using just enough tension to conform to the shape of the leg, ending just below the knee (figure 3).
- (4) Cut or hand-tear excess material and tape to keep the bandage in place (figure 4).


Layer 2 - Compression Layer:

- (1) Begin wrapping Layer 2 on the top of the foot, at the base of the fifth digit in a spiral technique (figure 5).
- (2) Hold the bandage with thumb and stretch until the ovals become circles. Begin spiral wrapping across the top of the foot, then proceed to wrap two figure of 8's around the back of the heel/ankle, completely covering the heel and foot. Continue to maintain even tension by stretching the visual indicators from ovals to circles (figure 6).

- (3) Continue to spiral wrap up the leg, covering Layer 1 with 50% overlap while maintaining even tension (figure 7).
- (4) Cut or hand-tear excess material (figure 8).
- (5) Gently press over entire application to ensure a strong bond & reduce slippage (figure 9).
- (6) Apply the nylon stocking over the system for ease of movement and patient comfort in footwear, under clothes and on bed sheets (figure 10).

Removal/Disposal

- CoFlex TLC Zinc (*Standard & Lite*) may be removed by unwrapping or with bandage scissors.
- Dispose of product per your facilities protocol.

Storage

For best results, store CoFlex TLC Zinc (*Standard & Lite*) in its original packaging at room temperature until use. Avoid excessive heat and humidity. Refer to the label on each box for expiration date.

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