

OVIK Health

COFLEX® UMAFD

Universal Medic Absorbent Foam Dressing



COFLEX® UMAFD is a durable trauma bandage designed for situations requiring an absorbent dressing, strength for compression, and high heat/water resistance. It is a non-latex, non-sterile, self-adhering bandage with an attached absorbent foam pad. The cohesive bandage allows the wrap to stick to itself but not to other material or skin with EasyTear technology intended for quick and easy application. This dressing is designed to be worn for up to three days.

Indications for Use:

CoFlex UMAFD is intended for use as an all-in-one trauma/combat bandage to compress and support with quick, easy application in high-stress situations.

- Great alternative to existing combat bandages
- Quick and easy to apply in high stress situations
- High absorption and retention of fluids (equivalent to 29 gauze pads)
- Flexible, controlled compression
- Able to withstand high temperatures (heat tested at 130° > 90 days)
- Made in the USA

Ideal for:

- Military - medics, soldiers
- First responders
- Emergency medical services
- ER / trauma centers
- Commercial and industrial sites
- First-aid kits

COFLEX UMAFD	Size	Package
9340UM-018	4" x 3 yd.	18 rolls/case
9340UM-100	4" x 3 yd.	100 rolls/case

Available in black only

