

OVIK Health™

COFLEX®

Flexible Cohesive Bandages

**PORTFOLIO
OVERVIEW**



Ovik
HEALTH



Improving Patient Lives. Together.

OVIK Health's CoFlex® cohesive bandages help patients and caregivers when a self-adhering bandage is desired to secure, protect, or compress. It is a gentle solution whenever adhesive is not wanted.

From proprietary latex-free formulation to EasyTear technology, OVIK Health provides the right solution for every clinical task.



QUALITY



VARIETY +
CUSTOMIZATION



COHESIVE
EXPERTISE



UNMATCHED
INNOVATION



EASYTEAR
TECHNOLOGY



CLINICAL
SUPPORT



BREATHABLE +
LIGHTWEIGHT

THE COFLEX DIFFERENCE

OVIK Health's CoFlex cohesive bandages are designed to simplify and streamline the healthcare clinicians job.

CoFlex SELECTION GUIDE

IMPROVES patient comfort with lightweight and breathable material

SECURES primary dressings and devices on patients

PROTECTS patients and clinicians with latex-free options

DELIVERS controlled compression that does not constrict over time

- Easy to apply
- Stays in place
- EasyTear technology
- Compresses to control swelling
- Conforms to area of application

LATEX

COFLEX

- Economical
- Variety of colors
- Quick stick
- Patented technology



NON_LATEX

COFLEX NL

- Clean hand-tear
- Variety of colors
- Prints
- Quick stick
- Available sterile
- Patented technology



COFLEX MED

- Clean hand-tear
- Variety of colors
- Quick stick
- Available sterile
- Water-resistant
- Patented technology



COFLEX LF2

- Clean hand-tear
- Variety of colors
- Quick stick
- Available sterile
- Conformable foam layer
- Water-resistant
- Patented technology





**BEST
PROTECTION**



COFLEX LF²

OVIK's CoFlex LF² is ideal for high-risk patient populations and surgical patients. CoFlex LF² is designed to provide comfort to patients with sensitive skin, eliminate the risk associated with latex sensitivity and protect the sterile field in surgical procedures.

- Soft and comfortable – gentle on the skin
- Ideal for OR application
- Water resistant for arthroscopic procedures**
- Sticks directly to a wet primary dressing
- EasyTear technology*
- Controlled compression – will not constrict or loosen over time
- Stays in place – does not shift or unravel with use



CoFlex LF² is available in color pack or tan!

Color Pack (CP)

Red, Blue, Neon Pink,
Neon Green, Purple,
Light Blue



Tan (TN)



Not made with natural rubber latex



**BEST
PERFORMANCE**



COFLEX NL

OVIK's CoFlex NL is a worry-free choice designed for use by everyone, everywhere. CoFlex NL provides excellent self adherence, excellent controlled compression, and excellent versatility to meet the multiple needs of an entire healthcare facility.

- Largest latex-free product offering with more than 90 options
- EasyTear technology
- Controlled compression - will not constrict or loosen over time
- Stays in place - does not shift or unravel with use
- Breathable and lightweight
- Made in the USA



CoFlex NL is available in a variety of colors and prints!

Color Pack (CP)

Red, Blue, Neon Pink, Neon Green, Purple, Light Blue



Rainbow Pack (RB)

Teal, Blue, White, Purple, Red, Green



Kids Pack (KP)

Bumble Bees, Hearts, Cars & Trucks, Dinosaurs - 2", 3", 4"



Smiley Pack (SC)

Smiley on Light Blue, Lavender, Yellow, Neon Pink - 2" only



Tan (TN)



Blue (BL)



Neon Green (NG)
2" only

Red (RD)

3" only



Light Blue (LB)
4" only



Not made with natural rubber latex



**BEST
VALUE**



COFLEX & COFLEX MED

OVIK Health's CoFlex & CoFlex Med offer economical solutions when price is paramount. CoFlex and CoFlex Med contain latex and provide exceptional self adherence while providing controlled compression.

- Controlled compression – will not constrict or loosen over time
- Stays in place – does not shift or unravel with use
- Versatile and can be used for a wide variety of applications
- Breathable and lightweight
- EasyTear technology*
- Water resistant for arthroscopic procedures**



CoFlex and CoFlex Med are available in a variety of color packs and tan!

CoFlex Options:

Color Pack (CP)

Red, Blue, Neon Pink,
Neon Green, Purple,
Light Blue



Rainbow Pack (RB)

Teal, Blue, White,
Purple, Red, Green



Tan (TN)



Blue (BL)



Red (RD)
3" only

CoFlex Med Options:

Color Pack (CP)

Red, Blue, Neon Pink,
Neon Green, Purple,
Light Blue



Tan (TN)



Not made with natural rubber latex




HCPCS CODES

(Healthcare Common Procedure Coding System Codes)


Product	Size	HCPCS Code
CoFlex CoFlex Med CoFlex NL CoFlex LF ²	1" x 5 yds. 1.5" x 5 yds. 2" x 5 yds.	A6453
CoFlex CoFlex Med CoFlex NL CoFlex LF ²	3" x 5 yds. 4" x 5 yds.	A6454
CoFlex CoFlex Med CoFlex NL CoFlex LF ²	6" x 5 yds.	A6455

*Centers for Medicare & Medicaid Services (CMS) has verified that the HCPCS codes listed above have been assigned to Andover Healthcare, Inc.
(Updated July 2015)


APPLICATION AND REMOVAL TIPS



1. The CoFlex family of bandages will not loosen or tighten over time, therefore you should wrap the bandage at desired tension and compression.




2. Upon completion of wrapping, grasp bandage between fingers and remove stretch - hand-tear* or cut the bandage to desired length. Relax the last 2" of bandage and then press down.



3. Seal the end of the bandage wrap with your finger, lightly meshing the top of the bandage into the bottom layer. This will provide long staying power.

Removal: CoFlex bandages may be removed by unwrapping or with bandage scissors.

YOU CONTROL THE TENSION



1. No to low tension
2. Medium tension
3. High tension

ORDERING INFORMATION

Catalog #	Size	Rolls/Case	Layer 2 (Stretched)		
CoFlex LF² †6 inch available in Tan only					
9100XX-030	1" x 5 yds.	30	9300(S)-024 <small>sterile/tan</small>	3" x 5 yds.	24
9100(S)-024 <small>sterile/tan</small>	1" x 5 yds.	24	9400XX-018	4" x 5 yds.	18
9150XX-048	1.5" x 5 yds.	48	9400(S)-018 <small>sterile/tan</small>	4" x 5 yds.	18
9200XX-036	2" x 5 yds.	36	9400L-018 <small>tan</small>	4" x 6.5 yds.	18
9200(S)-024 <small>sterile/tan</small>	2" x 5 yds.	24	9600XX-012†	6" x 5 yds.	12
9300XX-024	3" x 5 yds.	24	9600(S)-012 <small>sterile/tan</small>	6" x 5 yds.	12
CoFlex NL †6 inch available in Tan only, ‡1.6 yard available in Tan only					
5100XX-030	1" x 5 yds.	30	5600XX-012†	6" x 5 yds.	12
5150XX-048	1.5" x 5 yds.	48	5600(S)-012 <small>sterile/tan</small>	6" x 5 yds.	12
5200XX-036	2" x 5 yds.	36	5100X1-030‡	1" x 1.6 yds.	30
5300XX-024	3" x 5 yds.	24	5150X1-048‡	1.5" x 1.6 yds.	48
5300(S)-024 <small>sterile/tan</small>	3" x 5 yds.	24	5200X1-036‡	2" x 1.6 yds.	36
5400XX-018	4" x 5 yds.	18	5300X1-024‡	2" x 1.6 yds.	24
5400(S)-018 <small>sterile/tan</small>	4" x 5 yds.	18	5400X1-018‡	4" x 1.6 yds.	18
5400L-018 <small>tan</small>	4" x 6.5 yds.	18			
CoFlex †6 inch available in Tan only					
3100XX-030	1" x 5 yds.	30	3400XX-018	4" x 5 yds.	18
3150XX-048	1.5" x 5 yds.	48	3400L-018 <small>tan</small>	4" x 6.5 yds.	18
3200XX-036	2" x 5 yds.	36	3600XX-012†	6" x 5 yds.	12
3300XX-024	3" x 5 yds.	24			
CoFlex Med †6 inch available in Tan only					
7100XX-030	1" x 5 yds.	30	7400XX-018	4" x 5 yds.	18
7150XX-048	1.5" x 5 yds.	48	7400(S)-018 <small>sterile/tan</small>	4" x 5 yds.	18
7200XX-036	2" x 5 yds.	36	7400L-018 <small>tan</small>	4" x 6.5 yds.	18
7200(S)-036 <small>sterile/tan</small>	2" x 5 yds.	36	7600XX-012†	6" x 5 yds.	12
7300XX-024	3" x 5 yds.	24	7600(S)-012 <small>sterile/tan</small>	6" x 5 yds.	12
7300(S)-024 <small>sterile/tan</small>	3" x 5 yds.	24			

When ordering, please indicate catalog # and color code. For example, 9100TN-030 (XX = color code).

ovikhealth.com | 800.432.6686

©2023 OVIK Health, LLC. All rights reserved.

*EasyTear applicable for CoFlex Med, CoFlex NL and CoFlex LF². Scissors required for CoFlex.

**Water-resistant applicable for CoFlex Med and CoFlex LF².

OVIKCohesivePO r0 2023-10

