


OVIK Health™

SPORTS MEDICINE

PORTFOLIO OVERVIEW



Ovik
HEALTH

THE POWERFLEX DIFFERENCE

For the Athlete

- Supports muscles and joints
- Prevents injury
- Provides post-injury stability

For the Application

- Eliminates need for prewrap and adhesive spray
- Increases speed of application
- Compatible as the foundation for taping systems

For the Athletic Trainer


- Streamlines ordering by reducing materials needed from 3 to 1
- Allows for modification to address each athlete's needs
- Delivers consistent application

OVIK Health's cohesive wraps and tapes are designed to empower athletes to perform at their best while simplifying and streamlining the athletic trainers job.

In addition to these products, we support athletic trainers with literature and videos, sharing expertise on taping methods for specific injury types.

Backed by an organization with more than 150 years of experience in textile design, chemical science and materials manufacturing, we strive to improve clinical outcomes and enhance athletic performance while also helping our customers to increase efficiency, reduce waste and optimize their budgets.

Visit our website for more information or to request product samples.



WHICH TAPE OR WRAP SHOULD I CHOOSE?

BEST Overall



ULTIMATE TAPING SYSTEM

Combining the strength and support of **POWERFLEX + POWERTAPE**

Are YOU using PowerFlex as a base?

- 23 lbs. tensile strength
- Sweat & water resistant
- Replaces the need for prewrap and adhesive spray

Combine it with PowerTape for ultimate support.

- Versatile adhesive tape
- Less than 1% stretch
- Cotton free - won't absorb sweat or water

BEST Strength

POWERFLEX & POWERTAPE

(cohesive athletic tape)

BEST Budget

POWERFLEX & VICTORYTAPE

(cohesive athletic tape)

BEST Speed

POWERSPEED & POWERTAPE

(30% stronger than PowerFlex - apply less)

BEST Reinforcement

POWERFAST

(Apply over any system for additional support without restricting movement)

Build-Your-Own Taping System:

BASE LAYER (straight on skin)

POWERFLEX Self-Adhering Sports Wrap

POWERSPEED The Ultimate Speed Tape

POWERFAST Flexible Adhesive Stretch

TOP LAYER (tape)

POWERTAPE Cohesive Athletic Tape

VICTORYTAPE Cohesive Athletic Tape

SELF-ADHERING SPORTS WRAP

POWERFLEX

POWERFLEX and POWERFLEX LF

PowerFlex® empowers stronger athletes with its superior performance, wide variety of applications and its compatibility as the foundation of taping systems. This strong, fabric-based wrap tears easily by hand and sticks to itself for quick wrapping and longer staying power.

Use PowerFlex as a base to streamline ordering. Replacing prewrap, adhesive spray, and heel & lace pads, it reduces materials needed from 3 to 1!

- 23 lbs. tensile strength
- Sweat & water resistant
- Eliminates need for prewrap and adhesive spray
- Supports muscles & joints
- Provides injury prevention and post-injury stability for inversions, eversions and high ankle sprains
- Will not slip, slide or migrate
- #1 wrap in professional football & college sports



| PowerFlex Catalog # | Size | Rolls/Case | Notes |
|---------------------|----------------|------------|--------------------------------|
| 3710XX-048 | 1" x 6 yds. | 48 | Available in non-latex (White) |
| 3715XX-032 | 1.5" x 6 yds. | 32 | |
| 3720XX-024 | 2" x 6 yds. | 24 | |
| 3725XX-016 | 2.75" x 6 yds. | 16 | |
| 4725XX-048 | 2.75" x 6 yds. | 48 | |
| 4730XX-048 | 3" x 6 yds. | 48 | |
| 3730XX-016 | 3" x 6 yds. | 16 | |
| 3740XX-012 | 4" x 6 yds. | 12 | |
| 3760XX-008 | 6" x 6 yds. | 8 | |

Color Options:



*Yellow and Tan available in 2" only.

When ordering, please indicate catalog # and color code. For example, 3710WH-030 (XX = color code).

Additional

STRETCH LAYER OPTIONS

POWERSPEED POWERFAST

POWERSPEED

PowerSpeed® is a high strength cohesive athletic wrap designed for speed taping or a stronger spat. It delivers maximum support for a modified taping system that won't slip/migrate from sweat or water. With easy hand-tear, PowerSpeed has quick application.

- 30% stronger than PowerFlex
- Supports sprains & strains



POWERFAST

PowerFAST® is a strong, flexible adhesive stretch tape made with synthetic fabric. Unlike other adhesive stretch tapes, it won't absorb water so it remains thin & comfortable, rather than getting bulky. PowerFast also retains its elasticity for lasting support, plus its strong adhesive won't pull hair.

- Synthetic construction - won't absorb water
- Thin and comfortable - no bulk



| PowerSpeed Catalog # | Size | Rolls/Case |
|----------------------|----------------|------------|
| PS3720XX-024 | 2" x 6 yds. | 24 |
| PS3725XX-016 | 2.75" x 6 yds. | 16 |
| PS3730XX-016 | 3" x 6 yds. | 16 |

| PowerFAST Catalog # | Size | Rolls/Case |
|---------------------|-------------|------------|
| 4220XX-024 | 2" x 5 yds. | 24 |
| 4230XX-016 | 3" x 5 yds. | 16 |

Color Options:



When ordering, please indicate catalog # and color code. For example, PS3720WH-024 (XX = color code).

Additional
**TAPE LAYER
OPTIONS**

POWERTAPE VICTORYTAPE

POWERTAPE

PowerTape® is a cohesive athletic tape that sticks only to itself or to PowerFlex, creating one of the strongest combination taping systems in the field.

- 41 lbs. tensile strength
- Less than 1% stretch



VICTORYTAPE

VictoryTape® is an easier-to-tear and economical version of PowerTape to apply over PowerFlex for a long-lasting taping system.

- 38 lbs. tensile strength
- Sweat & water resistant



MOLESKIN

Moleskin is a tan-colored flannel with a napped cotton backing.



| PowerTape Catalog # | Size | Rolls/Case | VictoryTape Catalog # | Size | Rolls/Case |
|---------------------|----------------|------------|-----------------------|----------------|------------|
| ACP130XX-1" | 1" x 15 yds. | 48 | ACP120XX-1" | 1" x 15 yds. | 48 |
| ACP130XX-1.5" | 1.5" x 15 yds. | 32 | ACP120XX-1.5" | 1.5" x 15 yds. | 32 |
| ACP130XX-2" | 2" x 15 yds. | 24 | ACP120XX-2" | 2" x 15 yds. | 24 |

Color Options:



White (WH) Black (BK)

When ordering, please indicate catalog # and color code. For example, ACP130WH-1" (XX = color code).

| Moleskin Catalog # | Size | Rolls/Case | Colors |
|--------------------|--------------|------------|----------|
| 040-02-016* | 2" x 25 yds. | 16 | tan only |
| 040-03-016* | 3" x 25 yds. | 16 | tan only |

*Rolls are individually packaged

POWERFLEX AFD

PowerFlex® AFD is an all-in-one foam pad dressing and cohesive bandage that helps get players back in the game. It provides fast, effective treatment for wounds, lacerations and turf burns.

- Single use product
- Individually packaged
- Perfect for game time kits



| PowerFlex AFD Catalog # | Size | Rolls/ Case |
|-------------------------|---------------|-------------|
| 3722XX-016 | 2" x 2.5 yds. | 16 |
| 3744XX-008 | 4" x 2.5 yds. | 8 |
| 3744SB-008 (sterile)* | 4" x 2.5 yds. | 8 |

Color Options:



When ordering, please indicate catalog # and color code. For example, 3722WH-016 (XX = color code).

*Sterile available in black only.

Additional resources for our sports medicine options including product information sheets, application videos and case studies are available at www.ovikhealth.com.

